

Reactive Red M8B
C.I. Reactive Red 11

| | | |
|--------------------|---|--------|
| Light | : | 4 |
| Washing | : | 4-5 |
| Perspiration | : | 2 |
| Substantivity | : | High |
| Reactivity | : | High |
| Dischargeability | : | Poor |
| Hypochlorite | : | 4W |
| Solubility at 20°C | : | 70 g/l |

1% Cotton Dyeing
4% Cotton Dyeing
1% Hosiery Dyeing
4% Hosiery Dyeing



Reactive Magenta MB
C.I. Reactive Violet 13

| | | |
|--------------------|---|--------|
| Light | : | 4-5 |
| Washing | : | 3 |
| Perspiration | : | 4-5 |
| Substantivity | : | High |
| Reactivity | : | High |
| Dischargeability | : | Fair |
| Hypochlorite | : | 1 |
| Solubility at 20°C | : | 70 g/l |

1% Cotton Dyeing
4% Cotton Dyeing
1% Hosiery Dyeing
4% Hosiery Dyeing

